

Discover Anti-Aging and Skin Care with the Scalar Wave Laser including Main Unit and Probes



INTRODUCING: The Scalar Wave Laser (SWL)

Facial & Cosmetic
Regeneration and beautification.

A close-up photograph of a woman's face, focusing on her eyes and nose. Her skin appears slightly red and irritated, possibly due to acne or inflammation.

Anti-Aging
Slow the effects of time.

A close-up photograph of a woman's face, focusing on her forehead and eye area. She has visible wrinkles and a slightly concerned expression, suggesting the effects of aging.

Experience Cold Laser Facial Rejuvenation, Scar Reduction, Anti-Aging & Relief from Pain & Inflammation all in one laser.

Reduce or Eliminate Wrinkles and Acne with Ease. The Scalar Wave Laser works as an ideal skin care systems, which offers superior power, performance and surface area coverage with multiple frequencies. Utilizing 3 optimal wavelengths, a digitally blended light is achieved, which delivers superior results. Cold Laser Therapy (CLT) promotes the production of collagen and elastin, firms and tones skin, boosts immunity, eases inflammation, reduces redness from blemishes, abrasions and acne, heals and softens scars and allows for a natural “face lift” without drugs, injections or surgery.

NOW YOU CAN EXPERIENCE Pain Relief, Facial Rejuvenation & Anti-Aging with CLT

All natural, pain free and totally safe. Cold laser therapy and rejuvenation devices have been used actively in the medical, beauty and anti-aging field around the world for over 25 years. This stuff really works! Practitioners tell us that this is the best photo-facial and cold laser healing equipment they have ever used and they can't believe the prices compared to other professional units they have purchased, which are not as convenient, effective and easy to use, or approved by the FDA for consumer use.

Rejuvenate your body from the inside out. The Scalar Wave Laser (SWL) is also an ideal anti-aging, pain management and all around healing laser. It has 11 pre-programmed frequency setting that support the rejuvenation of skin, organs, glands, muscles, bones, tendons, ligament, meridians and charkas by working on a cellular level to donate photonic energy, while promoting anti-inflammatory and analgesic effects. The SWL gently and effectively unwinds stress, tension and old holding patterns out of the body and clears cellular memory to support perfect healing of the whole body.

TESTIMONIALS:

* I am getting younger! I heard about the unwinding process and began to do it daily. Within about 6 weeks people began to say, "have you had work done?" (cosmetic surgery) The work was on my thymus (anti-aging gland) appears to have turned back the clock!

– Cathy S.

* I am an RN who twisted my knee while walking my dog. After several months of trying many things to relieve the pain, I had no relief. Doctors told me that knee replacement was my only option. It was then that I attended a Quantumwave Gathering and had the opportunity to use the laser. During the gathering the pain was so improved that I was able to go dancing that evening. I had a couple of additional sessions the following month. It has been 9 months and the knee pain has not returned.

– Adele A

* The scars I have from chiggers (six years ago) are almost gone from using the scalar wave laser.

- Dawn R.

* We have a lady in Manitoba that burnt herself so badly, she was scheduled for skin grafts and extreme reconstructive surgery. After 2 weeks of 3 hours/day she is back to her job as a music teacher, and the doctors are completely dumbfounded and amazed. Her hands were burnt to the bone and 2 weeks later, you can hardly tell aside from some pink skin.

- Jalquoi A.

* I have been using the laser for about a month. I have lost weight and the swelling in my left leg has gone down and circulation increased. I have more energy and awareness. Also experienced pain relief in my hips. I feel very relaxed after I use the laser. Thank you!

– Irving G.

* 1.5 years ago I had an appointment for a facelift but cancelled at the last minute. Something didn't feel right inside. After purchasing the laser I began seeing a remarkable difference in my face. A girlfriend I hadn't seen in several months told me "You look 10 years younger. What are you doing?" She promptly purchased a laser also!

– Karen W.

Read Many More Testimonials Online:

<http://www.discoverlasers.com/scalar-wave-laser-testimonials>

The Scalar Wave Laser has OTC over the counter FDA Clearance



All natural, pain free and totally safe

Cold laser therapy and rejuvenation devices have been used actively in the medical, beauty and anti-aging field around the world for over 25 years. This stuff really works! Practitioners tell us that this is the best photo-facial and cold laser healing equipment they have ever used and they can't believe the prices compared to other units they have purchased, which are not as convenient, effective and easy to use.

Rejuvenate your body from the inside out and the outside in

The Scalar Wave Laser (SWL) is also an ideal anti-aging, pain management and healing laser. It has 11 pre-set frequency settings that support the rejuvenation of skin, organs, glands, muscles, bones, tendons, ligament, meridians and charkas by working on a cellular level to donate photonic energy, while promoting anti-inflammatory and analgesic effects. The SWL gently and effectively unwinds stress, tension and old holding patterns out of the body and clears cellular memory to support perfect healing. It combines the most advanced cold laser technology with state of the art quantum scalar wave technology to address all facets of health, dis-ease and addictions.

While most lasers are limited to square wave-forms and analog technology, **the Scalar Wave Laser is fully digital and utilizes multiple wave forms including:**

- Sine waves
- Scalar waves and
- Square waves

CONTENTS

Scalar Wave Laser Overview	1-3
Scalar Wave Laser Testimonials	2
Table of Contents	4
Introduction	5
The Main Scalar Wave Laser (SWL) System	5
Treatment Time for a full facial takes between 24 and 30 minutes	6
Contra-Indications	6
Rejuvenate Your Skin for Life - Circulation and Collagen	6
Treatment Instructions for the Scalar Wave Laser	7
Assembly	7
Operation of the Scalar Wave Laser & Choosing a Frequency Setting	8
The 11 SWL Essentials Settings Explained	8
Operation of Your Scalar Wave Laser with Pulsar & Viowave Probes	11
Customizing Features & How to View Remaining Battery Power	11
General Laser Light Therapy Instructions	12
Using Skin Care Products with the Scalar Wave Laser	12
Anti-Aging Laser Light Therapy Treatment Instructions	12
Care and Maintenance of Your SWL	14
To reduce the risk of injury or electrical shock	14
Precautions	14
Risks and Benefits	15
Scalar Wave Laser can be used on animals as well as humans	16-17
Scalar Wave Business Opportunity	18
Disclaimer and Copyright	18

Introduction

Welcome to the new frontier of Laser Light Therapy. Now you can experience advanced healing and anti-aging with cold laser or soft laser therapy utilizing the Scalar Wave Laser. You'll be happy to know that these devices are engineered for many years of trouble free use with only minimum maintenance and care. These lasers are FDA Cleared with Over The Counter (OTC) Clearance, which means that they are extremely safe for regular use and that anyone can buy & use them without a prescription. These lasers are manufactured in the USA to medical quality standards.

The Main Scalar Wave Laser (SWL) System

The main SWL is an ultra-portable, hand-held, rechargeable, cordless, multi-wavelength, resonating cluster laser with 3 optional stimulating, pulsed laser probes. Use of the probes will be discussed later in this manual. The SWL combines the most advanced cold laser technology with state of the art quantum scalar technology to significantly support optimum health and banish disease on a cellular level.

While most lasers are limited to square wave-forms and analog technology, **the Scalar Wave Laser is fully digital and utilizes multiple wave forms to effectively saturate the cells with photonic energy at multiple depths of penetration.**

With the SWL you get much more for less money. A total of 160mW of power spread over 2.5 inches surface area with 16 5mW true laser diodes and 20 4mW super violet LEDs. The SWL contains a winning combination of 8 visible red 650 nm and 8 invisible infrared 780 nm laser diodes, plus 20 violet 405 nm LEDs. With the SWL turned on, if you were to look indirectly at the aperture or laser face (don't look without glasses or blinking rapidly), you would see 8 laser diodes shining red, 20 LEDs shining violet and another 8 laser diodes that appear colorless or invisible. Do not worry, the colorless diodes are working, you simply cannot see the light they emit because infrared wavelengths are invisible to the human eye.

Laser Light Therapy is remarkably safe and effective on the skin, with no reported negative side effects, except occasional dryness of the skin. Follow treatment guidelines and do not over treat, because this can result in dry skin, in which case simply reduce the number of treatments per week. **Treatments of one to three times per week are usually sufficient.**

Treatment time for a full facial takes between 24 and 30 minutes.

While performing the treatment, the face should be clean with all makeup removed. Most users divide the face into six sections with a vertical line dividing the face from top to bottom and two horizontal lines going from left to right through the tip of the nose and the other just over the eyes. Then each of the 6 sections is treated for a total of four (4) minutes. Lightly apply the SWL to the skin to minimize light from escaping out the sides. You may gently and slowly move the laser in the section (called scanning), or you may keep it stationary for some minutes and then move it to another adjoining area.

Contra-Indications

Contra-Indications are conditions for which medical supervision or consultation is suggested before administering treatment. Using Light for facials or body treatments is an emerging technology and the interactions between the light and certain drugs and maladies have not been fully studied. Nor is it advisable to use directly over the fetus of a pregnant woman or directly over an active cancer, even though some people report skin cancer diminishing.

Certain medicines and products can cause light sensitive skin.

Under these conditions a medical doctor's permission and supervision is recommended before performing treatments. Contra-indications include:

- Epilepsy (pulsing red light can initiate seizures)
- Drugs that increase photosensitivity
- Cosmetic injections such as Botox – wait two weeks before beginning light therapy

Rejuvenate Your Skin for Life – Circulation and Collagen

Now you can easily support the production of collagen and elastin, which firms and tones skin, heals scars and essentially allows for a natural “face lift” without drugs, injections or surgery. The SWL uses the ultimate combination

of red, infrared and violet to support healing and rejuvenation of the skin cells and tissues. The SWL is in a class of medical devices using optimum wavelengths and frequencies to effect positive changes in the human body and on the skin.

Generally speaking, these changes are in cellular activity, blood flow, circulation, reduced inflammation and the production of collagen and elastin beneath the skin. Microcirculation and vasodilation are improved, swelling and acne is reduced and healing is accelerated, even in the case of old scars or after surgical procedures, with wound healing and when dealing with herpes simplex and zoster. The visible changes result in healthier, younger looking skin. A diminished appearance of wrinkles and fine lines; smaller pores; and improved color. Endorphins, enzymes and pain blockers are released to reduce pain and accelerate healing.

Friends will notice the change, even though they may not be able to say exactly what that difference is. The SWL treatment combines well with other treatments such as Microderm-abrasion, chemical peels and plastic surgery. It enhances the benefits gained from other skin care products and treatments, healing is much faster, and inflammation is markedly reduced, making other treatments more effective. Even make-up goes on more easily and often less is needed.

Treatment Instructions for the Scalar Wave Laser

Assembly

The SWL comes standard with an internal rechargeable lithium battery, an international power supply that can be used around the world and a US plug for ease of use in plugging into any outlet in North America.

USE ONLY THE POWER SUPPLY THAT COMES WITH THE SWL.

To operate, simply plug the power supply into an AC outlet using the correct plug or power cord as required for your country. Plug the line with the connector to the bottom of the SWL and then fully charge for 2-4 hours until all four bars show on the battery life user screen under “Settings”. The laser is often received fully charged. You can check it by turned on the laser (push the STAR button) and then scroll down to settings and then scroll to the right once. You will see time and bars. If there are less than 4 bars, then charge it until all 4 bars are shown. If you ever find yourself in a situation where the battery has run out of power and you wish to continue using the laser, simply plug the SWL in with the included power cord, and you can use it while it is charging.

Operation of the Scalar Wave Laser & Choosing a Frequency Setting

The SWL has multiple modes of operation, with 11 preset settings to choose from. For your convenience the settings have been labeled with names for easy identification.

A couple of them that you may wish to try are as follows: Unwind, Soothe, Vital and Wave. Simply turn on your SWL by pushing down on the STAR button located in the center of the star image on the hand piece. Once turned on the LCD screen will light up with a welcome screen and then quickly advance to the starting screen where you can choose from: Launch, Design and Settings. To Launch, simply scroll with the arrow to the right once and you will see Essentials and scroll one more time to the right and you will see the selection of the 11 Essential Settings or Modes. Now you can simply scroll to the right once more to activate the default setting, which is Unwind, or you can scroll down to your choice of Setting and then activate it by scrolling to the right. It will turn on and stay on for 4 minutes unless you purposely change the time interval setting under Settings.

The 11 SWL Essentials Settings Explained:

Essentials represents the 11 core settings used in the SWL. Each setting is a proprietary combination of frequencies that set the SWL apart from other lasers. They were specially chosen by the inventor, after extreme research and testing. They basically represent protocols for unwinding, pain relief, detoxification, regeneration, pleomorphic, immune boosting, organ, gland, emotion and chakra balancing, lymph and quantum activation applications. All frequency settings will transfer through for use with the probes as well. Following is a summary.

The Unwind Setting is the most commonly used program or mode because it unwinds cell memory, stress and tension from the body, glands, organs, Chakras, cells and energy systems. It is great to use on old injuries or scars and on any situation, which has stress, tension and contraction. Any time there is sickness, injury, stress or trauma, the cells and tissues go into contraction. In order for healing to take place and for balance to be restored, the cells and tissues need to be brought back to a state of neutrality, which is what unwind does best. Using Unwind for everything is ok. The basic gist of this program is the dissolving of contracted cellular memory and the return to the natural state of ease, which is neutrality. This is

the easiest program to start with and consequently, the default program that can be used initially for any application, because it works on the principle that all injury or dis-ease causes contraction. To dissolve contraction and bring cells and tissues back to a state of neutrality, supports and accelerates the healing process. If you are lazy or confused about which program to use, just use Unwind and either place the laser over the area of concern or follow the unwinding protocol to support the entire body, bringing ease, gentle energy and balance into all of the systems so that they can work efficiently to rebuild and rejuvenate.

The Quantum Setting is about shifting energy back into the quantum state. It works in a similar fashion to the Unwinding Program only it does so on the energetic and emotional level. It is a very popular program that works great on all cells, organs, Chakras and anywhere that is tense and contracted or in need of emotional clearing. Reported to be good for lower back pain among others. The Relieve Program is great for pain relief and injuries anywhere on the body and also helps balance the neurotransmitter sequencing with the brain to change the pain signaling into something different such as natural opiate signaling. Also good for probe use on ear points and body pains. Contains the popular Serotonin frequencies, which are also good for addictions, weight loss & dieting.

The Pleo Setting is a detox and immune activation program. It is based on research of the pleomorphic somatidian life cycle – and in particular the regenerative state of the somatid. Our goal here is to unwind, normalize and activate the somatid to its regenerative state. Pleo is based on higher dimensional quantum DNA physics and the principle of subtle regenerative coherence. This program is great for injuries, pain and any type of cellular regeneration including the support of hair growth and thickening. It's also a good one for using with probes over acupoints and for facial rejuvenation. I particularly like to use it 2-3 inches below the naval for 20 min. to cleanse the entire blood stream and over the liver to support detoxification.

The Rejuv Setting is designed to activate the lymph system. We have twice as much lymph in our body as we do blood and no heart/pump to move it, so lymph stagnation is common, especially in those who do not get regular exercise. Rejuv helps to induce the electrons into a high spin state. When electrons shift into a high spin state they phase shift into an ether state such as lymph, cranial sacral, brain, sexual and glandular substances. This program is a higher dimensional wellness mode to induce this high spin alchemical state.

The Soothe Setting is a very soothing program for injuries, pain and stress. It works great by dissolving tension. *It's also good for use with the probes over acupoints on the body, face and ears.* Good for doing facial rejuvenation, anti-wrinkle and acne treatment processes as it thoroughly supports the anti-inflammatory process.

Removes the inflammation and creates a supportive healing environment for the skin and old injuries and even scars. The Unity Program is ideal for neutralizing anything and everything, supporting a shift into unity. It works wonderfully through the violet laser light in the main unit and the violet probe for processing in-form-ation into the Quantum State. Works nicely with the Violet probe to scan the meridians and clear the energy in the Chakras and open the crown Chakra to the abundance of the spirit and a direct channel to in-form-ation that we can receive from higher consciousness.

The Vital Setting is about activating the inner cellular energy. It is great for activating the immune response, activating oxygen and electrons, detoxification and clarifying energy, regeneration and wellness. I like to use this one and Pleo any time I feel the need for an extra boost to kick off a cold, flu or sore throat. I also like this one for bringing vitality into my organs such as Pancreas, Liver, Lungs, Heart, Spleen, Intestines and Kidneys. Contains the popular Ozone frequency.

The Cohere Setting is ideal for bringing coherence to the body, Chakras (energy centers) and all of its systems. The body is designed in and of coherence and all health and or dis-ease is simply the reflection of how coherent we are as energy beings body/field and how efficiently we can process in-form-ation. Contains the popular Chakra frequencies.

The Wave Setting is for dissolving all cells and phenomena into one wave essence and then allowing that wave to expand into consciousness. Great for working around the head, and top 3 Chakra energy centers. Supports the flowering of consciousness and opening to the divine. Also supportive of increased hair growth. Good for supporting the frontal lobe and anterior singulet of the brain through use on the forehead and sphenoid and also good for use over the tip of the nose. These are all valuable points that support better health and balance. Extensive brain research at prestigious universities in the US have shown that regular stimulation of the frontal lobe & anterior singulet (through sphenoid) brain sections (10 min per day), may help to lower blood pressure, decrease pain, diminish depression, reduce stress and lower the risk of Alzheimer's, Dementia & mental decline.

To Launch simply arrow twice to the right and then scroll down through the ESSENTIALS and choose from the 11 Settings: Unwind, Quantum, Relieve, Pleo, Activate, Rejuv, Soothe, Unity, Vital, Cohere and Wave. Once you have scrolled down to the setting that you desire, then arrow to the right once to activate the laser. The default time setting is 4 minutes after which (if the sound is turned on) a soft tone will sound to let you know to move to the next area. At this time, in order to continue for an additional time with the same setting you simply push the arrow button to the right once more and you can resume with the same setting you were previously using. To pause the laser in the middle of a treatment simply arrow once to the left and then to resume again arrow once to the right.

Operation of Your Scalar Wave Laser with Pulsar & Viopulsar Probes

If using the SWL with a Pulsar Probe attached, the main SWL base unit becomes the power supply and frequency generator for the probe. The Red 650 nm, 100 mW probe and the Violet 405 nm, 100 mW probe are both popular anti-aging lasers for supporting skin healing, smoothing wrinkles, natural face lift, treatment of acupoints, & easing acne & inflammation. Probes are ideal for delivering a more concentrated healing energy faster, to desired points of concern, which require increased energy and therapeutic advantage. The Violet probe is good for anti-aging & supports the length of the Telomerase enzymes & activates DNA sequencing. Use probes for 1-3 minutes in each region or directly over lines, wrinkles, acne, blisters, sores, herpes, burns, dermatitis, psoriasis & more. It will quickly support the cells in reestablishing youthfulness & vitality.

Customizing Features & How to View Remaining Battery Power

If you wish to increase or decrease the treatment time interval between shut off, turn the sound on or off, adjust the contrast or backlighting or *view the amount of battery power remaining in the your laser*, simply return to the starting screen and arrow down to Settings. At Settings, arrow once to the right and then arrow up or down to choose the field you wish to adjust. When adjusting the time interval bear in mind that the time will return to the 4 minute default once the laser has been shut off.

To change Session length, arrow to the right once, and then the “EDITING

Settings” screen will appear. At this point use the key pad to enter the desired number of minutes, and then arrow back to the left to save the new entry and then arrow back up to “Launch” and proceed as previously instructed.

General Laser Light Therapy Instructions

A course of anti-aging facial laser light therapy involves approximately six weeks with a treatment one to three times per week.

Following the series of initial treatments, maintenance treatments will be required every 30 to 60 days to maintain the benefits, or more frequently for potentially greater results. The actual time between maintenance treatments will vary depending on your age and skin. You will have to determine the frequency of maintenance treatments from your own experience. **Do not over treat. It could dry out the skin.**

Using Skin Care Products with the Scalar Wave Laser

When performing anti-aging light therapy treatments, it is highly advisable to use a serum of hyaluronic acid. This highly active and powerfully hydrating serum is extremely complimentary in laser light therapy, as it drives it deeply into the cells for enhanced rejuvenation. Aging skin is linked to intracellular dehydration. Combining laser light therapy with high quality natural skin care products will enhance the properties of the product. You may safely use laser light therapy with skin products that are otherwise safe to use on the skin. Make certain that the product emulsifies into the skin and does not leave a barrier that prevents the laser light from penetrating. **You may safely use skin care products and sun screens after laser light therapy if desired.**

Anti-Aging Laser Light Therapy Treatment Instructions

Before you begin, remember to remove all makeup and jewelry from the areas to be treated. Cleanse these areas thoroughly.

Next, divide your face into 6 Sections, splitting the face into left and right halves and then into 3 horizontal areas. Imagine your face is a clock with 6 sections: 1, 3, 5, 7, 9 and 11 o'clock. You will be starting at 1 o'clock and after 4 minutes progressing to 3 o'clock and then to 5, 7, 9 and 11 o'clock.

Begin with Section 1: Choose your Setting such as Soothe or Unwind, and then place your SWL gently against your skin. You may leave it in one place or move your SWL in small circles, slowly but continuously over the designated area. You may also use a red or violet probe in place of or in addition to the main unit for additional benefits, especially for deeper lines.

A tone will sound after 4 minutes, signaling you to move on to the next section.

Move to section 3.

When treating under the eye, gently follow the cheekbone. Avoid shining laser light into your eyes by keeping them closed or by covering them with a slice of cooling cucumber.

Again you will treat Section 3 for 4 minutes, before hearing the signal tone.

Move to section 5.

Continue treating sections 7, 9, and 11 as previously instructed.

Total treatment will take 24 minutes.

Press the STAR button on the back of the SWL until the unit turns off. You will know it is turned off when the LCD screen is no longer on.

Do not block the treatment areas with heavy creams or lotions.

You may use any moisturizer or sun block after completing your full SWL treatment.

Care and Maintenance of Your SWL

The SWL has been built for service and reliability, and should give many years of service if a few service tips are observed.

Turn your SWL off when not in use. We also recommend turning it off between treatments in order to allow the diodes to cool and in order to save on battery power.

Clean outside of hand piece with an alcohol swab or damp (not wet) cloth after each use.

Do not use dishwasher or autoclave to clean or sterilize. Do not put in hot or boiling water.

Do not place hard plastic case near open flame or other heat source.

CAUTION! To reduce the risk of injury or electrical shock:

Avoid using the SWL while standing in water (including bath and shower).

Avoid dropping the device in water. Water can short-circuit the electrical parts of your SWL, which could result in shock or damage.

Do not attempt to use the SWL if it has fallen in water or other liquids. Water may damage the parts resulting in electric shock.

Do not use the SWL if the power supply has a frayed cord or the prongs are bent or loose.

Do not disassemble the SWL as this may cause damage, malfunction, or fire, which could result in electrical shock. There are no user-serviceable parts inside the unit. Disassembly will void the warranty.

CAUTION! Precautions:

To help with treatments and get the best results proceed as follows:

Always clean your SWL after each use. Unplug the SWL, use a damp (not wet) cloth or isopropyl alcohol swab, and wipe the surfaces.

To avoid damage to your SWL and keep your warranty in force:

Use only the power supply provided with the SWL device. Other power units could result in damage to the SWL.

Never place the power cord near heated surfaces such as stoves or curling irons. The cord may melt which could lead to a short circuit or loss of power to your SWL.

Never scratch, mark, or otherwise damage the light emitting surface or lenses. This could affect or lower the light output or lead to electric shock.

Do not place in water or attempt to sterilize with heat. It could cause the unit to malfunction or fail to operate.

Risks and Benefits

The SWL device will reduce the signs of aging if you follow the directions given in this Instruction Manual. For persons with lines, wrinkles, or enlarged pores, Quantum light therapy offers a method for improving the appearance of the skin without medication, or injections.

You could experience, however, one or more of the following conditions: Mild skin redness or pain in the area of treatment. Less than 1% of users have reported either. If either of these occurs and the effect lasts longer than you feel comfortable with or if the redness or pain lasts longer than 24 hours, contact your doctor. The red and infrared and violet technology used in the SWL is well within the safe limits as established by International standards.



The Scalar Wave Laser can be used on animals as well as humans

They are popular for use on horses, dogs and cats. The main benefits include:

- * Reduces inflammation**
- * Relieves pain**
- * Regenerates tissue**
- * Promotes circulation**

What sets this cold laser apart from any other laser is its ability to unwind and clear cell memory. As with all therapies, veterinary permission should always be sought before treating any animal. There are hundreds of research papers from academic sources to back the efficacy of cold laser therapy for use with animals.

Cold Laser Therapy has been used with competition horse training and recovery for over 30 years. It works well for race-horses, show jumpers, show horses and ponies, endurance horses, carriage horses, polo ponies, quarter horses and old retired horses and ponies. By using the SW Laser on wounds and muscle tissue, it speeds up the body's own natural healing process several times over. This helps to prevent unsightly scar tissue and the white hairs you so often see on old scar sites. The treatment only takes a short time, can be done frequently and has no known side effects, allowing your horse, dog or cat to return to health quickly.

Low level lasers have proven to be beneficial and aid the following conditions in horses, dogs, cats & most animals.

- * **General.** Stress/distress, muscle pain, trauma, swellings, joint pain & arthritis, stiffness.
- * **Wounds.** Deep cuts, proud flesh, scars, swellings, bruising, edema, lumps.
- * **Operation sites.** Over stitches, staples, pins, pain relief, faster wound healing, less scarring.
- * **Back.** Pain, cold back, sore back, arthritis, bruises, strains, fibrotic myopathies, muscle pain, kissing spine.
- * **Hip trouble.** Sacroiliac symptom.
- * **Skin.** Acute infections, chronic infections, eczema, sunburn, cracked heels, mud fever, rain scold, haematomas, non-healing ulcers.
- * **Legs.** Tendon injuries, tendonitis, bowed tendons, ligament injuries, inflammation, adhesions, torn fibres, strains, splints, arthritis, Spavins, ringbone, bruising, sore shins.
- * **Feet.** Navicular, cracks & splits, laminitis, bruised soles, over-reach wounds, abscess, thrush.
- * **Head.** Sinuses, laryngeal paralysis, congested tear ducts, tooth abscess, sore gums & cheeks.



Horses, Dogs and Cats all Love the Scalar Wave Laser. They relax, breath more fully and respond quickly with pain relief, more energy, happiness and playfulness.

**Come and join the Cold Laser Revolution.
Call Laser Specialist Kalon Prensky for support or to order
additional lasers at preferred customer prices: 1-808-870-0348**

The Scalar Wave Laser is revolutionizing the way people and practitioners approach healing, skin care, pain management, rehabilitation, anti-aging, addictions, smoking cessation, weight control, animal care and much, much more. Now you can get on board with this rapidly growing light therapy revolution of complimentary health care. Now FDA Cleared!

Disclaimer and Copyright:

The Scalar Wave Laser System is FDA indicated for temporary relief of minor muscle and joint pain, arthritis and muscle spasm, relieving stiffness, promoting relaxation of muscle tissue, and temporarily increasing local blood circulation.

The information and protocols provided within this report are for informational purposes only. Neither the manufacturer nor distributors of the devices discussed, nor the author make any claims for use on people to diagnose, treat or cure any disease or disorder. Laser products, testimonials and information are not to be considered a replacement for any treatment, therapy or medication prescribed by a healthcare practitioner or Doctor. No health claims are intended or implied. Statements made are for informational purposes only and have not been evaluated or approved by the U.S. Food & Drug Administration.

If you have a medical condition, please see a qualified physician. The manufacturers and sales representatives of the lasers shown and described make no human claims regarding their products except for those that it has been FDA cleared for. Kalon Prensky, independent distributor assumes no liability for any of the SWL products sold. In the case of product defect or warranty, etc., customers must contact manufacturer directly. All information contained within is protected by copyright law. No part of it may be copied or reproduced without written consent by Kalon Prensky. Distribution of it in any form is illegal.

© 2006-2015 by Kalon Prensky